

SA100+HF

MINXRAY[®]


High Frequency Radiology System
For Small or Mixed Animal Practice



Big dog? No problem.



Horse skull in the field?
No problem either!



Compact and versatile system to meet a wide range of imaging needs

MinXray's unique **SA100+HF** system provides high functionality for in-clinic and ambulatory small and large animal imaging. This system offers powerful, state-of-the-art x-ray technology packaged with an easy to use space-saving table.

SA100+HF

- requires just 10 square feet of floor space
- includes the **XRT300** table with cassette tray and 6:1 ratio, 103 lpi stationary grid
- installs easily and operates on incoming voltage of 110 - 240 VAC 50/60 Hz
- provides easy mounting and removal of portable x-ray device
- is fully supported by MinXray's superior 5-year limited warranty
- is backed by an unequalled loaner program providing the end user with years of worry free ownership

